



OUR PRODUCTS

OUR CLIENTS

ABOUT US

CONTACT US

SUPPORT

ARTICLES

NEWS & OFFERS

PRODUCT SPECIFICATIONS

TERMS & CONDITIONS



The Active Birth Pool

- The Active Birth Pool's innovative ergonomic design provides unparalleled comfort, support and safety for mothers and midwives
- The pool's unique organic shape is psychologically approachable and aesthetically pleasing, helping mothers to feel more relaxed in the hospital environment.
- A rim level safety seat allows emergency evacuation and in-pool monitoring
- The pool's free standing design allows midwives access to mothers from all sides in comfort
- Our skirting panel makes the pool more attractive, faster to clean and easier to maintain
- A Wide rim enables mothers and midwives to support themselves comfortably
- Our Bespoke fast flow drainage system empties the pool in less than five minutes
- Handrails provide support for squatting and other upright positions
- An Internal Seat provides support during labour and birth, and can act as a step assisting entrance and exit from the pool
- The pool's gently textured non-slip surface provides secure stable footing

[Learn More >](#)

